

## **CORNING® 1L POLYCARBONATE ERLLENMEYER FLASK WITH DIP TUBE, WITH 0.2 MM PTFE FILTER, MLL**

Corning Erlenmeyer Flasks are now available as closed system solutions. They are sterile and arrive individually double-bagged and preassembled. The preassembled transfer caps have two ports. One port ends in a 0.2 µm filter, and the other port is chemically resistant, heat sealable flexible tubing, ending in a male Luer lock. The dip tube reaches all the way to the bottom of the flask for easy aseptic transfer of your liquid handling processes.

- Made from optically clear polycarbonate
- Sterilized by gamma radiation and nonpyrogenic

**Cat. #:** Corning 11440

**Qty./Pk:** 1 / Pk

**Qty./Cs:** 4 / Cs

**Brand:** Corning®

**Cap Material:** Polypropylene

**Cap Color:** Orange

**Sterile:** Yes

**Feature:** Disposable

**Capacity:** 1000 mL

**Flask Style:** Erlenmeyer

**Top Style:** Threaded

**Graduation Interval:** 50 mL



**Note:** This product was over-purchased from Corning. It is 70% off Corning list price.