

## **COSTAR® 2 ML ASPIRATING PIPETS, POLYSTYRENE, WITHOUT GRADUATIONS, INDIVIDUALLY WRAPPED, STERILE, 1/BAG, 1000/CASE**

### **Costar® 2 mL Aspirating Pipets, Polystyrene, Without Graduations, Individually Wrapped, Sterile, 1/Bag, 1000/Case**

Costar® aspirating pipets are sterile, without graduations and unplugged for aseptically aspirating liquid by using vacuum suction.

- One piece polystyrene construction
- Sterile and nonpyrogenic

**Cat. #:** Corning 9186

Qty./Pk: 100 / Pk

Qty./Cs: 1000 / Cs

Brand: Corning®

Capacity: 2 mL

Pipet Style: Aspirating

Tip Shape: Round

Feature: Disposable

Packaging: Individually wrapped

Length: 277 mm (approx.)



**Note:** This product was over-stocked from Corning. It is 50% off Corning list price.