

CYNOMOLGUS ADAM9 PROTEIN, HIS TAG**Cat.#:** 12118**Product Name:** Cynomolgus ADAM9 Protein**Size:** 10 µg, 50 µg and 100 µg**Synonyms:** CORD9;MCMP;MDC9;Mltng**Target:** ADAM9**UNIPROT ID:** A0A2K5X4X8**Description:** Recombinant Cynomolgus ADAM9 protein with C-terminal 6xHis tag**Background:** This gene encodes a member of the ADAM (a disintegrin and metalloprotease domain) family. Members of this family are membrane-anchored proteins structurally related to snake venom disintegrins, and have been implicated in a variety of biological processes involving cell-cell and cell-matrix interactions, including fertilization, muscle development, and neurogenesis. The protein encoded by this gene interacts with SH3 domain-containing proteins, binds mitotic arrest deficient 2 beta protein, and is also involved in TPA-induced ectodomain shedding of membrane-anchored heparin-binding EGF-like growth factor. Several alternatively spliced transcript variants have been identified for this gene. [provided by RefSeq, Jul 2010]**Species/Host:** HEK293**Molecular Weight:** The protein has a predicted molecular mass of 74.9 kDa after removal of the signal peptide. The apparent molecular mass of cADAM9-His is approximately 55-100 kDa due to glycosylation.**Molecular Characterization:** ADAM9(Ala29-Gly698) 6×His tag**Purity:** The purity of the protein is greater than 85% as determined by SDS-PAGE and Coomassie blue staining.**Formulation & Reconstitution:** Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.**Storage & Shipping:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.




Figure 1. Cynomolgus ADAM9 Protein, His Tag on SDS-PAGE under reducing condition.