

TNFSF11 MAB (DENOSUMAB)

Cat.#: 28066

Product Name: Anti-TNFSF11 (Denosumab) Monoclonal Antibody

Synonyms: CD254;hRANKL2;ODF;OPGL;OPTB2;RANKL;sOdf;TNLG6B;TRANCE

Description: Anti-TNFSF11 (denosumab biosimilar) Monoclonal Antibody

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Applications: Flow Cyt

Recommended Dilutions: Flow Cyt 1:100

Host Species: Homo sapiens

Isotype: IgG2

Purification: Purified from cell culture supernatant by affinity chromatography

Species Reactivity: Human TNFSF11

Constituents: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

Storage & Shipping: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. HEK293 cell line transfected with irrelevant protein (red histogram) and human TNFSF11 protein (blue histogram) were surface stained with Anti-TNFSF11 (denosumab biosimilar) mAb $1\mu\text{g/ml}$ followed by Alexa 488-conjugated Goat anti-Human IgG secondary antibody.