

CGRPR MAB (ERENUMAB)

Cat.#: 28157

Product Name: Anti-CGRPR(Erenumab) Monoclonal Antibody

Synonyms: CALCRL

Description: Anti-CGRPR(Erenumab) Monoclonal Antibody

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Applications: Flow Cyt

Recommended Dilutions: Flow Cyt 1:100

Host Species: Homo sapiens

Isotype: IgG2

Purification: Purified from cell culture supernatant by affinity chromatography

Species Reactivity: Human CGRPR

Constituents: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

Storage & Shipping: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. Flow cytometry analysis with $1\mu\text{g}/\text{mL}$ Anti-CGRPR(Erenumab biosimilar) mAb (28157) on Expi293 cells transfected with Human CGRPR&RAMPI protein (Blue histogram) and Expi293 transfected with Human CGRPR protein (Red histogram)) or Human RAMPI protein(Orange histogram).



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
