

FALCON® 96-WELL CLEAR FLAT BOTTOM NOT TREATED CELL CULTURE MICROPLATE, WITH LID, INDIVIDUALLY WRAPPED, STERILE, 50/CASE

Falcon® 96-well Clear Flat Bottom Not Treated Cell Culture Microplate, with Lid, Individually Wrapped, Sterile, 50/Case

- Crystal-grade virgin polystyrene
- Nonpyrogenic
- Non-treated polystyrene
- Convenient, peel-open medical-style packaging for individual product presentation
- Individual trays are PET (Code 1) and recyclable.
- Falcon 96-well microplates are compatible with automation systems and meet standard ANSI/SBS footprint dimensions for microplates
- Working volume of Media is 0.1-0.2 mL

Cat. #: Falcon® 351172

Qty./Pk: 1 / Pk

Qty./Cs: 50 / Cs

Brand: Falcon®

Plate Format: 96-well

Plate Color: Clear

Well Bottom: Flat

Well Volume: 0.37 mL

Cell Growth Area: 0.32 cm²


Recommended Working Volume: 40-275 µL

Surface Treatment: Not Treated

Sterile: Yes

Lids Included: Yes

Storage: Ambient storage conditions



Note: This product was over-stocked from Corning. It is 50% off Corning list price.