

## **FALCON® 96-WELL CLEAR FLAT BOTTOM NOT TREATED CELL CULTURE MICROPLATE, WITH LID, INDIVIDUALLY WRAPPED, STERILE, 50/CASE**

### **Falcon® 96-well Clear Flat Bottom Not Treated Cell Culture Microplate, with Lid, Individually Wrapped, Sterile, 50/Case**

- Crystal-grade virgin polystyrene
- Nonpyrogenic
- Non-treated polystyrene
- Convenient, peel-open medical-style packaging for individual product presentation
- Individual trays are PET (Code 1) and recyclable.
- Falcon 96-well microplates are compatible with automation systems and meet standard ANSI/SBS footprint dimensions for microplates
- Working volume of Media is 0.1-0.2 mL

**Cat. #:** Falcon® 351172

**Qty./Pk:** 1 / Pk

**Qty./Cs:** 50 / Cs

**Brand:** Falcon®

**Plate Format:** 96-well

**Plate Color:** Clear

**Well Bottom:** Flat

**Well Volume:** 0.37 mL

**Cell Growth Area:** 0.32 cm<sup>2</sup>

**Recommended Working Volume:** 40-275 µL

**Surface Treatment:** Not Treated

**Sterile:** Yes

**Lids Included:** Yes

**Storage:** Ambient storage conditions



**Note:** This product was over-stocked from Corning. It is 50% off Corning list price.