

FALCON® 14 ML ROUND BOTTOM HIGH CLARITY PP TEST TUBE, GRADUATED, WITH SNAP CAP, STERILE, 25/PACK, 500/CASE**Falcon® 14 mL Round Bottom High Clarity PP Test Tube, Graduated, with Snap Cap, Sterile, 25/Pack, 500/Case**

Falcon® 17x100mm, 14 mL high-clarity polypropylene (PP) round bottom test tubes have a 3,000 RCF rating. Each tube features printed blue graduations and a white writing patch. The tubes come with polyethylene dual-position snap caps, which offer both vented and fully closed positions. The sterile tubes are sold 25 to a pack and 500 to a case. RNase/DNase-free.

- Widely used and referenced in laboratory protocols
- Dual position snap-cap prevents sample loss with ease of handling
- Sterile, nonpyrogenicity tested to less than 0.1 EU/ mL

Cat. #: Falcon® 352059

Qty./Pk: 25 / Pk

Qty./Cs: 500 / Cs

Brand: Falcon®

Tube Style: Test

Cap Style: Snap cap

Cap Color: Natural

Graduations: Yes

Capacity: 14 mL

Sterile: Yes

O.D.: 17 mm

Product Material: Polypropylene

Application: Diagnostics



Note: This product was over-stocked from Corning. It is 50% off Corning list price.