

## **FALCON® 48-WELL CLEAR FLAT BOTTOM TC-TREATED CELL CULTURE PLATE, WITH LID, INDIVIDUALLY WRAPPED, STERILE, 50/CASE**

### **Falcon® 48-well Clear Flat Bottom TC-treated Cell Culture Plate, with Lid, Individually Wrapped, Sterile, 50/Case**

Flat bottom, with low evaporation lid. Individually wrapped, sterile, 50/case. Unique labyrinth lid, condensation rings, and deep-well design control contamination, reduce evaporation, and minimize edge effects. Reliable vacuum-gas plasma tissue culture treatment provides well-to-well and plate-to-plate consistency.

- Crystal-grade virgin polystyrene
- Nonpyrogenic
- Standard tissue culture (TC)-treated
- Convenient, peel-open medical-style packaging for individual product presentation
- Individual trays are PET (Code 1) and recyclable
- Working volume of Media is 0.5-0.8 mL
- Working volume of Trypsin is 0.05-0.08 mL (Your lab protocol may call for another proteolytic enzyme.)

**Cat. #: Falcon® 353078**

**Qty./Pk:** 1 / Pk

**Qty./Cs:** 50 / Cs

**Brand:** Falcon®

**Plate Format:** 48-well


**Well Shape:** Round

**Cell Growth Area:** 0.75 cm<sup>2</sup>

**Surface Treatment:** TC-Treated

**Sterile:** Yes

**Lids Included:** Yes



**Note:** This product was over-stocked from Corning. It is 50% off Corning list price.