

**FALCON® 2 ML SEROLOGICAL PIPET, POLYSTYRENE, 0.01 INCREMENTS,
INDIVIDUALLY PACKED, STERILE, 100/BOX, 1,000/CASE**

**Falcon® 2 mL Serological Pipet, Polystyrene, 0.01 Increments, Individually
Packed, Sterile, 100/Box, 1,000/Case**

For tissue culture usage or any non-toxic solution that requires aliquoted or aspirating/dispensing in a sterile environment. Nonpyrogenic pipet with 0.8 mL overflow capacity. Individually packed in a paper/plastic thermoformed package.

Cat. #: Falcon® 357507

Qty./Pk: 100 / Pk

Qty./Cs: 1000 / Cs

Brand: Falcon®

Capacity: 2 mL

Pipet Style: Serological

Feature: Disposable

Graduation Interval: 0.01 mL



Note: This product was over-stocked from Corning. It is 50% off Corning list price.