

**FALCON® 2 ML SEROLOGICAL PIPET, POLYSTYRENE, 0.01 INCREMENTS,  
INDIVIDUALLY PACKED, STERILE, 100/BOX, 1,000/CASE**

**Falcon® 2 mL Serological Pipet, Polystyrene, 0.01 Increments, Individually  
Packed, Sterile, 100/Box, 1,000/Case**

For tissue culture usage or any non-toxic solution that requires aliquoted or aspirating/dispensing in a sterile environment. Nonpyrogenic pipet with 0.8 mL overfill capacity. Individually packed in a paper/plastic thermoformed package.

**Cat. #:** Falcon® 357507

**Qty./Pk:** 100 / Pk

**Qty./Cs:** 1000 / Cs

**Brand:** Falcon®

**Capacity:** 2 mL

**Pipet Style:** Serological

**Feature:** Disposable

**Graduation Interval:** 0.01 mL



**Note:** This product was over-stocked from Corning. It is 50% off Corning list price.