

**FALCON® 2 ML ASPIRATING PIPET, POLYSTYRENE, WITHOUT GRADUATIONS,
INDIVIDUALLY WRAPPED, STERILE, 50/BAG, 200/CASE**

**Falcon® 2 mL Aspirating Pipet, Polystyrene, without Graduations,
Individually Wrapped, Sterile, 50/Bag, 200/Case**

Ideal for tissue culture usage or any non-toxic solution that requires aliquoted or aspirating/dispensing in a sterile environment. Nonpyrogenic, nonplugged pipet, individually wrapped in thermoformed paper/plastic wrap.

Cat. #: Corning 357558

Qty./Pk: 50 / Pk

Qty./Cs: 200 / Cs

Brand: Falcon®

Capacity: 2 mL

Pipet Style: Aspirating

Graduation Interval: 0.1 mL

Packaging: Bulk pack



Note: This product was over-stocked from Corning. It is 50% off Corning list price.