

**FALCON® 2 ML ASPIRATING PIPET, POLYSTYRENE, WITHOUT GRADUATIONS, INDIVIDUALLY WRAPPED, STERILE, 50/BAG, 200/CASE**

**Falcon® 2 mL Aspirating Pipet, Polystyrene, without Graduations, Individually Wrapped, Sterile, 50/Bag, 200/Case**

Ideal for tissue culture usage or any non-toxic solution that requires aliquoted or aspirating/dispensing in a sterile environment. Nonpyrogenic, nonplugged pipet, individually wrapped in thermoformed paper/plastic wrap.

**Cat. #:** Corning 357558

**Qty./Pk:** 50 / Pk

**Qty./Cs:** 200 / Cs

**Brand:** Falcon®

**Capacity:** 2 mL

**Pipet Style:** Aspirating

**Graduation Interval:** 0.1 mL

**Packaging:** Bulk pack



**Note:** This product was over-stocked from Corning. It is 50% off Corning list price.