

## **FISHERBRAND™ 96-WELL DEEPWELL™ POLYPROPYLENE MICROPLATES**

### **Fisherbrand™ 96-Well DeepWell™ Polypropylene Microplates**

Ideal for storage because of the low-binding characteristics of polypropylene

- Accommodate even larger working volumes; up to 1.0 mL or 2.0 mL storage volume per well
- Proteins and DNA will not adhere to PP, allowing for complete sample recovery.
- Natural color
- Nonsterile and sterile available

**Cat. #:** Fisherbrand™ 12-566-612

**Format:** Standard

**Color:** Natural

**Volume (Metric) Well:** 2000  $\mu$ L

**Volume (Metric) Working:** 1900  $\mu$ L

**Well Shape:** Round

**Array:** 8 x 12

**Material:** Polypropylene

**No. of Wells:** 96

**Sterility:** Sterile

**For Use With (Application):** For storage of low-binding characteristics of polypropylene

