

## **Product Description**

Pioneering GTPase and Oncogene Product Development since 2010

## FISHERBRAND™ 96-WELL DEEPWELL™ POLYPROPYLENE MICROPLATES

## Fisherbrand™ 96-Well DeepWell™ Polypropylene Microplates

Ideal for storage because of the low-binding characteristics of polypropylene

- Accommodate even larger working volumes; up to 1.0 mL or 2.0 mL storage volume per well
- Proteins and DNA will not adhere to PP, allowing for complete sample recovery.
- Natural color
- Nonsterile and sterile available

Cat.#: Fisherbrand™ 12-566-612; 12566612

Format: Standard

**Color:** Natural

Volume (Metric) Well: 2000 μL

Volume (Metric) Working: 1900 μL

Well Shape: Round

**Array:** 8 x 12

Material: Polypropylene

No. of Wells: 96 Sterility: Sterile

For Use With (Application): For storage of low-binding characteristics of

polypropylene



## **Product Description**

Pioneering GTPase and Oncogene Product Development since 2010

