

FISHERBRAND™ 96-WELL DEEPWELL™ POLYPROPYLENE MICROPLATES

Fisherbrand™ 96-Well DeepWell™ Polypropylene Microplates

Ideal for storage because of the low-binding characteristics of polypropylene

- Accommodate even larger working volumes; up to 1.0 mL or 2.0 mL storage volume per well
- Proteins and DNA will not adhere to PP, allowing for complete sample recovery.
- Natural color
- Nonsterile and sterile available

Cat. #: Fisherbrand™ 12-566-612

Format: Standard

Color: Natural

Volume (Metric) Well: 2000 µL

Volume (Metric) Working: 1900 µL

Well Shape: Round

Array: 8 x 12

Material: Polypropylene

No. of Wells: 96

Sterility: Sterile

For Use With (Application): For storage of low-binding characteristics of polypropylene

