

**HUMAN FZD10 PROTEIN, HFC TAG****Cat.#:** 11446**Product Name:** Human FZD10 Protein**Size:** 10 µg, 50 µg and 100 µg**Synonyms:** CD350;FZ-10;Fz10;FzE7;hFz10**Target:** FZD10**UNIPROT ID:** Q9ULW2**Description:** Recombinant Human FZD10 with C-terminal human Fc tag**Background:** This gene is a member of the frizzled gene family. Members of this family encode 7-transmembrane domain proteins that are receptors for the Wingless type MMTV integration site family of signaling proteins. Most frizzled receptors are coupled to the beta-catenin canonical signaling pathway. Using array analysis, expression of this intronless gene is significantly up-regulated in two cases of primary colon cancer. [provided by RefSeq, Jul 2008]**Species/Host:** HEK293**Molecular Weight:** The protein has a predicted molecular mass of 42.2 kDa after removal of the signal peptide. The apparent molecular mass of FZD10-hFc is approximately 35-55 kDa due to glycosylation.**Molecular Characterization:** FZD10(Ile21-Gly161) hFc(Glu99-Ala330)**Purity:** The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.**Formulation & Reconstitution:** Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.**Storage & Shipping:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.




Figure 1. Human FZD10 Protein, hFc Tag on SDS-PAGE under reducing condition.