

Product Description

Pioneering GTPase and Oncogene Product Development since 2010

HUMAN FZD4 PROTEIN, HFC TAG

Cat.#: 11414 Product Name: Human FZD4 Protein Size: 10 µg, 50 µg and 100 µg Synonyms: CD344;EVR1;FEVR;Fz-4;Fz4;FZD4S;FzE4;GPCR;hFz4 Target: FZD4 UNIPROT ID: Q9ULV1

Description: Recombinant Human FZD4 with C-terminal human Fc tag

Background: This gene is a member of the frizzled gene family. Members of this family encode seven-transmembrane domain proteins that are receptors for the Wingless type MMTV integration site family of signaling proteins. Most frizzled receptors are coupled to the beta-catenin canonical signaling pathway. This protein may play a role as a positive regulator of the Wingless type MMTV integration site signaling pathway. A transcript variant retaining intronic sequence and encoding a shorter isoform has been described, however, its expression is not supported by other experimental evidence. [provided by RefSeq, Jul 2008]

Species/Host: HEK293

Molecular Weight: The protein has a predicted molecular mass of 42.4 kDa after removal of the signal peptide. The apparent molecular mass of FZD4hFc is approximately 55-70 kDa due to glycosylation.

Molecular Characterization: FZD4(Phe37-Glu180) hFc(Glu99-Ala330) **Purity:** The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Formulation & Reconstitution: Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

Storage & Shipping: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010



Figure 1. Human FZD4 Protein, hFc Tag on SDS-PAGE under reducing condition.