

## HUMAN IL11 PROTEIN, HFC TAG

**Cat.#:** 11435

**Product Name:** Human IL11 Protein

**Size:** 10 µg, 50 µg and 100 µg

**Synonyms:** AGIF;IL-11

**Target:** IL11

**UNIPROT ID:** P20809

**Description:** Recombinant Human IL11 with C-terminal human Fc tag

**Background:** The protein encoded by this gene is a member of the gp130 family of cytokines. These cytokines drive the assembly of multisubunit receptor complexes, all of which contain at least one molecule of the transmembrane signaling receptor IL6ST (gp130). This cytokine is shown to stimulate the T-cell-dependent development of immunoglobulin-producing B cells. It is also found to support the proliferation of hematopoietic stem cells and megakaryocyte progenitor cells. Alternatively spliced transcript variants encoding distinct isoforms have been found for this gene. [provided by RefSeq, Jun 2012]

**Species/Host:** HEK293

**Molecular Weight:** The protein has a predicted molecular mass of 45.3 kDa after removal of the signal peptide. The apparent molecular mass of IL11-hFc is approximately 55–70 kDa due to glycosylation.

**Molecular Characterization:** IL11(Pro22–Leu199) hFc(Glu99–Ala330)

**Purity:** The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

**Formulation & Reconstitution:** Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

**Storage & Shipping:** Store at –20°C to –80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at –80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.



Figure 1. Human IL11 Protein, hFc Tag on SDS-PAGE under reducing condition.