

## HUMAN/MOUSE FGF-8B PROTEIN

**Cat.#:** 12080

**Product Name:** Human/Mouse FGF-8B Protein

**Size:** 10 µg, 50 µg and 100 µg

**Synonyms:** Fibroblast growth factor 8; Androgen-induced growth factor; Heparin-binding growth factor 8; AIGF; HBGF-8; FGF-8B

**Target:** FGF-8b

**UNIPROT ID:** P55075-3/P37237-2

**Description:** Recombinant Human/Mouse Fibroblast Growth Factor 8B is produced by our E.coli expression system and the target gene encoding Gln23-Arg215 is expressed.

**Background:** Fibroblast growth factor 8 (FGF-8) is a member of the fibroblast growth factor family. It is discovered as a growth factor essential for the androgen-dependent growth of mouse mammary carcinoma cells. Mouse FGF-8b shares 100% aa identity with human FGF-8b. FGF-8 is widely expressed during embryogenesis, and mediates epithelial-mesenchymal transitions. It plays an important role in the regulation of embryonic development, cell proliferation, cell differentiation and cell migration. It is required for normal brain, eye, ear, limb development during embryogenesis and normal development of the gonadotropin-releasing hormone (GnRH) neuronal system.

**Species/Host:** E.coli

**Molecular Weight:** 22.5 kDa

**Molecular Characterization:** Not available

**Purity:** Greater than 95% as determined by reducing SDS-PAGE.

**Formulation & Reconstitution:** Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

**Storage & Shipping:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.




Figure 1. Greater than 95% as determined by reducing SDS-PAGE.




Figure 2. Measured in a cell proliferation assay using BALB/c 3T3 cells. The ED50 for this effect is 21.87 ng/ml.