

HUMAN PGLYRP1 PROTEIN, HIS TAG**Cat.#:** 11576**Product Name:** Human PGLYRP1 Protein**Size:** 10 µg, 50 µg and 100 µg**Synonyms:** PGRP-S;PGRPS;PGLYRP;PGRP;TNFSF3L**Target:** PGLYRP1**UNIPROT ID:** O75594**Description:** Recombinant human PGLYRP1 protein with C-terminal 6xHis tag**Background:** Pattern receptor that binds to murein peptidoglycans (PGN) of Gram-positive bacteria. Has bactericidal activity towards Gram-positive bacteria. May kill Gram-positive bacteria by interfering with peptidoglycan biosynthesis. Binds also to Gram-negative bacteria, and has bacteriostatic activity towards Gram-negative bacteria. Plays a role in innate immunity.[UniProtKB/Swiss-Prot Function]**Species/Host:** HEK293**Molecular Weight:** The protein has a predicted molecular mass of 20.3 kDa after removal of the signal peptide. The apparent molecular mass of PGLYRP1-His is approximately 15–25 kDa due to glycosylation.**Molecular Characterization:** PGLYRP1(Gln22–Pro196) 6×His tag**Purity:** The purity of the protein is greater than 85% as determined by SDS-PAGE and Coomassie blue staining.**Formulation & Reconstitution:** Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.**Storage & Shipping:** Store at –20°C to –80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at –80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.




Figure 1. Human PGLYRP1 Protein, His Tag on SDS-PAGE under reducing condition.