

HUMAN TNFRSF1A PROTEIN, HFC TAG**Cat.#:** 11721**Product Name:** Human TNFRSF1A Protein**Size:** 10 µg, 50 µg and 100 µg**Synonyms:** Tumor necrosis factor receptor superfamily member 1A;NF-R1;NF-R1;NFR-I**Target:** TNFRSF1A**UNIPROT ID:** P19438**Description:** Recombinant human TNFRSF1A protein with C-terminal human Fc tag

Background: This gene encodes a member of the TNF receptor superfamily of proteins. The encoded receptor is found in membrane-bound and soluble forms that interact with membrane-bound and soluble forms, respectively, of its ligand, tumor necrosis factor alpha. Binding of membrane-bound tumor necrosis factor alpha to the membrane-bound receptor induces receptor trimerization and activation, which plays a role in cell survival, apoptosis, and inflammation. Proteolytic processing of the encoded receptor results in release of the soluble form of the receptor, which can interact with free tumor necrosis factor alpha to inhibit inflammation. Mutations in this gene underlie tumor necrosis factor receptor-associated periodic syndrome (TRAPS), characterized by fever, abdominal pain and other features. Mutations in this gene may also be associated with multiple sclerosis in human patients. [provided by RefSeq, Sep 2016]

Species/Host: HEK293

Molecular Weight: The protein has a predicted molecular mass of 46.6 kDa after removal of the signal peptide. The apparent molecular mass of TNFRSF1A-hFc is approximately 55-70 kDa due to glycosylation.

Molecular Characterization: TNFRSF1A(Leu30-Thr211) hFc(Glu99-Ala330)

Purity: The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Formulation & Reconstitution: Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

Storage & Shipping: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.




Figure 1. Human TNFRSF1A Protein, hFc Tag on SDS-PAGE under reducing condition.