

CD40 MAB (ISCALIMAB)

Cat.#: 28029

Product Name: Anti-CD40 (iscalimab) Monoclonal Antibody

Synonyms: CD40;Bp50;CDW40;MGC9013;TNFRSF5;p50

Description: Anti-CD40 (iscalimab) Monoclonal Antibody

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Applications: ELISA; Flow Cyt

Recommended Dilutions: ELISA 1:5000-10000; Flow Cyt 1:100

Host Species: Homo sapiens

Isotype: IgG1

Purification: Purified from cell culture supernatant by affinity chromatography

Species Reactivity: Human CD40

Constituents: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

Storage & Shipping: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).

Anti-CD40 (iscalimab biosimilar) mAb ELISA

0.2 µg of CD40, mFc-His Tagged protein per well




Figure 1. ELISA plate pre-coated by 2 µg/ml (100 µl/well) Human CD40, mFc-His tagged protein 11143 can bind Anti-CD40 Neutralizing antibody in a linear range of 0.64-80.0 ng/ml.




Figure 2. HEK293 cell line transfected with irrelevant protein (B) and human CD40 protein (A) were surface stained with anti-CD40 neutralizing antibody 1µg/ml (iscalimab) followed by Alexa 488-conjugated anti-human IgG secondary antibody.