

OX40 MAB (IVUXOLIMAB)

Cat.#: 28061

Product Name: Anti-OX40(Ivuxolimab)Monoclonal Antibody

Synonyms: TNFRSF4;OX40;CD134;OX40L receptor;ACT35;TXGP1L

Description: Anti-OX40(ivuxolimab biosimilar)Monoclonal Antibody

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Applications: Flow Cyt

Recommended Dilutions: Flow Cyt 1:100

Host Species: Homo sapiens

Isotype: IgG2

Purification: Purified from cell culture supernatant by affinity chromatography

Species Reactivity: Human OX40

Constituents: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

Storage & Shipping: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. Flow cytometry analysis with Anti-OX40 mAb 1 µg/ml (ivuxolimab) on Expi293 cells transfected with human OX40 (Blue histogram) or Expi293 transfected with irrelevant protein(Red histogram).

