

## IL7RA MAB (LUSVERTIKIMAB)

**Cat.#:** 28087

**Product Name:** Anti-IL7RA(Lusvertikimab) Monoclonal Antibody

**Synonyms:** CD127;CDW127;IL-7R-alpha;IL7RA;ILRA

**Description:** Anti-IL7RA(Lusvertikimab ) Monoclonal Antibody

**Background:** Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

**Applications:** Flow Cyt

**Recommended Dilutions:** Flow Cyt 1:100

**Host Species:** Humanized

**Isotype:** IgG4

**Purification:** Purified from cell culture supernatant by affinity chromatography

**Species Reactivity:** Human IL7RA

**Constituents:** Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

**Storage & Shipping:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. Flow cytometry analysis with 1 µg/mL Anti-IL7RA (lusvertikimab biosimilar) mAb (28087) on Expi293 cells transfected with Human IL7RA protein (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram).



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---