

IL7RA MAB (LUSVERTIKIMAB)

Cat.#: 28087

Product Name: Anti-IL7RA(Lusvertikimab) Monoclonal Antibody

Synonyms: CD127;CDW127;IL-7R-alpha;IL7RA;ILRA

Description: Anti-IL7RA(lusvertikimab) Monoclonal Antibody

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Applications: Flow Cyt

Recommended Dilutions: Flow Cyt 1:100

Host Species: Humanized

Isotype: IgG4

Purification: Purified from cell culture supernatant by affinity chromatography

Species Reactivity: Human IL7RA

Constituents: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

Storage & Shipping: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. Flow cytometry analysis with 1 µg/mL Anti-IL7RA (lusvertikimab biosimilar) mAb (28087) on Expi293 cells transfected with Human IL7RA protein (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram).

