

## NKG2A MAB (MONALIZUMAB)

**Cat.#:** 28088

**Product Name:** Anti-NKG2A(Monalizumab) Monoclonal Antibody

**Synonyms:** CD159A;NKG2;NKG2A

**Description:** Anti-NKG2A(monalizumab ) Monoclonal Antibody

**Background:** Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

**Applications:** Flow Cyt

**Recommended Dilutions:** Flow Cyt 1:100

**Host Species:** Humanized

**Isotype:** IgG4

**Purification:** Purified from cell culture supernatant by affinity chromatography

**Species Reactivity:** Human NKG2A

**Constituents:** Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

**Storage & Shipping:** Store at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$  for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at  $-80^{\circ}\text{C}$  (Avoid repeated freezing and thawing).




Figure 1. Flow cytometry analysis with 1  $\mu\text{g}/\text{mL}$  Anti-NKG2A (monalizumab biosimilar) mAb (28088) on Expi293 cells transfected with Human NKG2A protein and Human CD94 protein (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram).



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---