

Product Description

Pioneering GTPase and Oncogene Product Development since 2010

MOUSE TREM2 PROTEIN, HFC TAG

Cat.#: 12202 Product Name: Mouse TREM2 Protein Size: 10 µg, 50 µg and 100 µg Synonyms: TREM-2;Trem2a;Trem2b;Trem2c Target: TREM2

UNIPROT ID: Q99NH8

Description: Recombinant mouse TREM2 protein with C-terminal human Fc tag

Background: The protein encoded by this gene is part of the immunoglobulin and lectin-like superfamily and functions as part of the innate immune system. This gene forms part of a cluster of genes on mouse chromosome 17 thought to be involved in innate immunity. This protein associates with the adaptor protein Dap-12 and recruits several factors, such as kinases and phospholipase C-gamma, to form a receptor signaling complex that activates myeloid cells, including dendritic cells and microglia. In humans homozygous loss-of-function mutations in this gene cause Nasu-Hakola disease and mutations in this gene may be risk factors to the development of Alzheimer's disease. In mouse mutations of this gene serve as a pathophysiological model for polycystic lipomembranous osteodysplasia with sclerosing leukoencephalopathy (Nasu-Hakola disease) and for inflammatory bowel disease. Alternative splicing results in multiple transcript variants that encode different protein isoforms. [provided by RefSeq, Jan 2013]

Species/Host: HEK293

Molecular Weight: The protein has a predicted molecular mass of 42.9 kDa after removal of the signal peptide.

Molecular Characterization: Mouse TREM2(Leu19-Ser171) hFc(Glu99-Ala330)

Purity: The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Formulation & Reconstitution: Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

Storage & Shipping: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010



Figure 1. Mouse TREM2 Protein, hFc Tag on SDS-PAGE under reducing condition.