

## MASP2 MAB (NARSOPLIMAB)

**Cat.#:** 28135

**Product Name:** Anti-MASP2(Narsoplimab) Monoclonal Antibody

**Synonyms:** MASP-2

**Description:** Anti-MASP2(narsoplimab ) Monoclonal Antibody

**Background:** Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

**Applications:** Flow Cyt

**Recommended Dilutions:** Flow Cyt 1:100

**Host Species:** Homo sapiens

**Isotype:** IgG4

**Purification:** Purified from cell culture supernatant by affinity chromatography

**Species Reactivity:** Human MASP2

**Constituents:** Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

**Storage & Shipping:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. MASP2 protein is highly expressed on the surface of Jurkat cell membrane. Flow cytometry analysis with 15 µg/mL Anti-MASP2(narsoplimab biosimilar) mAb (28135) (Blue histogram) or isotype control mAb (Red histogram) on Jurkat cells.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---