

MASP2 MAB (NARSOPLIMAB)

Cat.#: 28135

Product Name: Anti-MASP2(Narsoplimab) Monoclonal Antibody

Synonyms: MASP-2

Description: Anti-MASP2(narsoplimab) Monoclonal Antibody

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Applications: Flow Cyt

Recommended Dilutions: Flow Cyt 1:100

Host Species: Homo sapiens

Isotype: IgG4

Purification: Purified from cell culture supernatant by affinity chromatography

Species Reactivity: Human MASP2

Constituents: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

Storage & Shipping: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. MASP2 protein is highly expressed on the surface of Jurkat cell membrane. Flow cytometry analysis with 15 µg/mL Anti-MASP2(narsoplimab biosimilar) mAb (28135) (Blue histogram) or isotype control mAb (Red histogram) on Jurkat cells.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
