

B7-H2 MAB (PREZALUMAB)

Cat.#: 28084

Product Name: Anti-B7-H2(Prezalumab) Monoclonal Antibody

Synonyms: ICOSLG;B7-H2;B7H2;B7RP-1;B7RP1;CD275;GL50;ICOS-L;ICOSL;LICOS;ICOS ligand

Description: Anti-B7-H2(prezalumab) Monoclonal Antibody

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Applications: Flow Cyt

Recommended Dilutions: Flow Cyt 1:100

Host Species: Homo sapiens

Isotype: IgG2

Purification: Purified from cell culture supernatant by affinity chromatography

Species Reactivity: Human B7-H2

Constituents: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

Storage & Shipping: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. Flow cytometry analysis with 15 µg/mL Anti-B7-H2 (prezalumab biosimilar) mAb (28084) on Expi293 cells transfected with Human B7H2 protein (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram).