

## B7-H2 MAB (PREZALUMAB)

**Cat.#:** 28084

**Product Name:** Anti-B7-H2(Prezalumab) Monoclonal Antibody

**Synonyms:** ICOSLG;B7-H2;B7H2;B7RP-1;B7RP1;CD275;GL50;ICOS-L;ICOSL;LICOS;ICOS ligand

**Description:** Anti-B7-H2(prezalumab ) Monoclonal Antibody

**Background:** Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

**Applications:** Flow Cyt

**Recommended Dilutions:** Flow Cyt 1:100

**Host Species:** Homo sapiens

**Isotype:** IgG2

**Purification:** Purified from cell culture supernatant by affinity chromatography

**Species Reactivity:** Human B7-H2

**Constituents:** Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

**Storage & Shipping:** Store at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$  for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at  $-80^{\circ}\text{C}$  (Avoid repeated freezing and thawing).




Figure 1. Flow cytometry analysis with 15  $\mu\text{g}/\text{mL}$  Anti-B7-H2 (prezalumab biosimilar) mAb (28084) on Expi293 cells transfected with Human B7H2 protein (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram).