

RAB18 PAB**Rab18 Polyclonal Antibody****Cat. #:** 21050**Gene Symbol:** RAB18**Description:** Anti-Rab18 Rabbit Polyclonal Antibody**Background:** Rab18 is a member of the Ras-related small GTPases. Rab18 plays a role in apical endocytosis/recycling.**Immunogen:** A synthetic peptide from the C-terminal end of human Rab18 protein**Applications:** ELISA, WB, IHC**Recommended Dilutions:**


ELISA: 1:100–1:5000

WB: 1:100–1:2000

IHC: 1:100–1:1000

Concentration: 0.2 mg/ml**Host Species:** Rabbit**Format:** Liquid**Clonality:** Polyclonal**Purity:** Purified from serum**Preservative:** no**Constituents:** PBS (without Mg^{2+} and Ca^{2+}), pH 7.4, 150 mM NaCl, 50% glycerol**Species Reactivity:** Rab18 antibody recognizes Rab18 proteins from vertebrates.**Storage Conditions:** Store at $-20^{\circ}C$. Avoid repeated freezing and thawing

Western blot:



WB: anti-Rab18 pAb

EC-1 (lane 1), HepG (lane 2) and HEK-293T (lane 3) whole-cell lysates were used for Western blots with anti Rab18 polyclonal antibody (cat # 21050).