

RAB18 PAB

Rab18 Polyclonal Antibody

Cat. #: 21050

Gene Symbol: RAB18

Description: Anti-Rab18 Rabbit Polyclonal Antibody

Background: Rab18 is a member of the Ras-related small GTPases. Rab18 plays a role in apical endocytosis/recycling.

Immunogen: A synthetic peptide from the C-terminal end of human Rab18 protein

Applications: ELISA, WB, IHC

Recommended Dilutions:

ELISA: 1:100–1:5000

WB: 1:100–1:2000

IHC: 1:100–1:1000

Concentration: 0.2 mg/ml

Host Species: Rabbit

Format: Liquid

Clonality: Polyclonal

Purity: Purified from serum


Preservative: no

Constituents: PBS (without Mg^{2+} and Ca^{2+}), pH 7.4, 150 mM NaCl, 50% glycerol

Species Reactivity: Rab18 antibody recognizes Rab18 proteins from vertebrates.

Storage Conditions: Store at $-20^{\circ}C$. Avoid repeated freezing and thawing

Western blot:



WB: anti-Rab18 pAb

EC-1 (lane 1), HepG (lane 2) and HEK-293T (lane 3) whole-cell lysates were used for Western blots with anti Rab18 polyclonal antibody (cat # 21050).