

## RAC1 MAB

### Rac1 Monoclonal Antibody

**Cat. #:** 26005

**Gene Symbol:** RAC1

**Description:** Anti-Rac1 Mouse Monoclonal Antibody

**Background:** Small GTPases are a super-family of cellular signaling regulators. Rac belongs to the Rho sub-family of GTPases that regulate cell motility, cell division, and gene transcription. In mammals, there are three Rac genes: Rac1, Rac2 and Rac3. GTP binding increases the activity of Rac, and the hydrolysis of GTP to GDP renders it inactive. GTP hydrolysis is aided by GTPase activating proteins (GAPs), while exchange of GDP for GTP is facilitated by guanine nucleotide exchange factors (GEFs).

**Immunogen:** Recombinant full length protein of Rac1

**Applications:** ELISA, WB, IHC

**Recommended Dilutions:** ELISA 1:1000-1:5000 WB 1:1000-1:2000

**Concentration:** 1 mg/ml

**Host Species:** Mouse

**Format:** Liquid

**Clonality:** Monoclonal

**Isotype:** IgG

**Purity:** Purified from ascites


**Preservative:** No

**Constituents:** PBS (without  $Mg^{2+}$  and  $Ca^{2+}$ ), pH 7.4, 150 mM NaCl, 50% glycerol

**Species Reactivity:** Rac1 antibody recognizes Rac1 from vertebrates.

**Storage Conditions:** Store at  $-20^{\circ}C$ . Avoid repeated freezing and thawing

### Western blot:



Western blot analysis of Rac1 proteins. MEF whole-cell lysate was used for western blots. The monoclonal Anti-Rac1 antibody was pre-incubated without (lane 1) or with (lane 2) the antigen protein (Rac1) to show the specificity