

GITR:TNFRSF18 MAB (RAGIFILIMAB)

Cat.#: 28072

Product Name: Anti-GITR:TNFRSF18(Ragifilimab) Monoclonal Antibody

Synonyms: AITR;GITR;TNFRSF18;CD357

Description: Anti-GITR:TNFRSF18(ragifilimab) Monoclonal Antibody

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Applications: Flow Cyt

Recommended Dilutions: Flow Cyt 1:100

Host Species: Humanized

Isotype: IgG4

Purification: Purified from cell culture supernatant by affinity chromatography

Species Reactivity: Human GITR

Constituents: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

Storage & Shipping: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. Flow cytometry analysis with Anti-GITR/TNFRSF18(ragifilimab) mAb 1 μ g/mL on Expi293 cells transfected with Human GITR (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram).