

SARS-COV-2 (2019-NCOV) S1 PROTEIN CTD, HFC TAG

Cat.#: 11300

Product Name: SARS-CoV-2 (2019-NCoV) S1 Protein CTD

Size: 10 µg, 50 µg and 100 µg

Synonyms: S1 protein CTD;pike protein S1 CTD;etaCoV S1-CTD;OVID-19

Target: S1 protein CTD

UNIPROT ID: P0DTC2

Description: Recombinant SARS-CoV-2 (2019-nCoV) S1 protein CTD with C-terminal human Fc tag

Background: SARS-CoV-2 (Severe Acute Respiratory Syndrome Coronavirus 2) also known as Covid19 (2019 Novel Coronavirus) is a virus that causes illnesses ranging from the common cold to severe diseases. The spike protein is a type I transmembrane protein containing two subunits, S1 and S2. S1 mainly contains a receptor binding domain (RBD), which accounts for recognizing the cell surface receptor, ACE2. S2 contains basic elements needed for the membrane fusion. Recent publications indicate that S1-RBD domain can induce virus neutralizing-antibody and T cell response.

Species/Host: HEK293

Molecular Weight: The protein has a predicted molecular mass of 47.9 kDa after removal of the signal peptide.

Molecular Characterization: S1 protein CTD(Asn334-Pro527) hFc(Glu99-Ala330)

Purity: The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Formulation & Reconstitution: Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

Storage & Shipping: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.




Figure 1. SARS-CoV-2 (2019-nCoV) S1 protein CTD, hFc Tag on SDS-PAGE under reducing condition.