

**SARS-COV-2 (KAPPA) NCOV S PROTEIN RBD PROTEIN, HFC TAG****Cat.#:** 11370**Product Name:** SARS-CoV-2 (Kappa) NCoV S Protein RBD Protein**Size:** 10 µg, 50 µg and 100 µg**Synonyms:** SARS-CoV-2 B.1.617/B.1.617.1/B.1.617.3 (Kappa) Spike RBD Protein**Target:** S protein RBD**UNIPROT ID:** P0DTC2**Description:** Recombinant SARS-CoV-2 RBD (L452R and E484Q) protein with C-terminal human Fc tag**Background:** SARS-CoV-2 (Severe Acute Respiratory Syndrome Coronavirus 2) also known as Covid19 (2019 Novel Coronavirus) is a virus that causes illnesses ranging from the common cold to severe diseases. The spike protein is a type I transmembrane protein containing two subunits, S1 and S2. S1 mainly contains a receptor binding domain (RBD), which accounts for recognizing the cell surface receptor, ACE2. S2 contains basic elements needed for the membrane fusion. Recent publications indicate that S1-RBD domain can induce virus neutralizing-antibody and T cell response.**Species/Host:** HEK293**Molecular Weight:** The protein has a predicted molecular mass of 51.2 kDa after removal of the signal peptide. The apparent molecular mass of RBD(L452R and E484Q)-hFc is approximately 55-70 kDa due to glycosylation.**Molecular Characterization:** S protein RBD(L452R and E484Q)(Arg319-Phe541) hFc(Glu99-Ala330)**Purity:** The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.**Formulation & Reconstitution:** Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.**Storage & Shipping:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.




Figure 1. SARS-CoV-2 (2019-nCoV) S protein RBD(L452R& E484Q), hFc Tag on SDS-PAGE under reducing condition.