

THERMO SCIENTIFIC™ NALGENE™ RAPID-FLOW™ STERILE FILTER STORAGE BOTTLES - 500 ML

Thermo Scientific™ Nalgene™ Rapid-Flow™ Sterile Filter Storage Bottles

Save time and minimize contamination with gamma-sterilized Thermo Scientific™ Nalgene™ Rapid-Flow™ Filter Storage Bottles. Use with Nalgene Bottle-Top filters or as replacement/extra receptacles for filter units. Excellent as single-use storage containers.

- Headspace allows for extra serum or other additives to be added to the filtrate
- Leakproof† screw caps to eliminate pH shift in bottle for longer life of stored media
- Sterilized by gamma radiation and individually bagged for a 5-year sterile shelf life
- Non-pyrogenic and non-cytotoxic

Cat. #: Thermo Scientific™ 4550500; Fisher Scientific 09-740-25C

Material: Virgon Clear Polystyrene Bottle, White Pigmented Virgin HDPE Closure

Model: Wide

No. per Case: 12

Volume (Metric): 500 mL

Sterility: Sterile

Mouth: Wide

No. per Pack: 12

Fit for Purpose: Sterile