

CXCR4 MAB (ULOCUPLUMAB)

Cat.#: 28101

Product Name: Anti-CXCR4(Ulocuplumab) Monoclonal Antibody

Synonyms: CD184;D2S201E;FB22;HM89;HSY3RR;LAP-3;LAP3;LCR1;LESTR;NPY3R;NPYR;NPYRL;NPYY3R;WHIM;WHIMS

Description: Anti-CXCR4(ulocuplumab) Monoclonal Antibody

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Applications: Flow Cyt

Recommended Dilutions: Flow Cyt 1:100

Host Species: Homo sapiens

Isotype: IgG4

Purification: Purified from cell culture supernatant by affinity chromatography

Species Reactivity: Human CXCR4

Constituents: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

Storage & Shipping: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. Flow cytometry analysis with 15 $\mu\text{g}/\text{mL}$ Anti-CXCR4 (ulocuplumab biosimilar) mAb (28101) on Expi293 cells transfected with Human CXCR4 protein (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram).