

## OSMR MAB (VIXARELIMAB)

**Cat.#:** 28140

**Product Name:** Anti-OSMR(Vixarelimab) Monoclonal Antibody

**Synonyms:** IL-31 receptor subunit beta;IL-31R subunit beta;IL-31R-beta;IL-31RB

**Description:** Anti-OSMR(vixarelimab ) Monoclonal Antibody

**Background:** Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

**Applications:** Flow Cyt

**Recommended Dilutions:** Flow Cyt 1:100

**Host Species:** Homo sapiens

**Isotype:** IgG4-lambda2

**Purification:** Purified from cell culture supernatant by affinity chromatography

**Species Reactivity:** Human OSMR

**Constituents:** Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

**Storage & Shipping:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. Flow cytometry analysis with 15 µg/mL Anti-OSMR (vixarelimab biosimilar) mAb (28140) on Expi293 cells transfected with Human OSMR protein (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram).