

OSMR MAB (VIXARELIMAB)

Cat.#: 28140

Product Name: Anti-OSMR(Vixarelimab) Monoclonal Antibody

Synonyms: IL-31 receptor subunit beta;IL-31R subunit beta;IL-31R-beta;IL-31RB

Description: Anti-OSMR(vixarelimab) Monoclonal Antibody

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Applications: Flow Cyt

Recommended Dilutions: Flow Cyt 1:100

Host Species: Homo sapiens

Isotype: IgG4-lambda2

Purification: Purified from cell culture supernatant by affinity chromatography

Species Reactivity: Human OSMR

Constituents: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

Storage & Shipping: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. Flow cytometry analysis with 15 $\mu\text{g}/\text{mL}$ Anti-OSMR (vixarelimab biosimilar) mAb (28140) on Expi293 cells transfected with Human OSMR protein (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram).