


ENPP3 MAB (XENCOR MAB)**Cat.#:** 28143**Product Name:** Anti-ENPP3(Xencor) Monoclonal Antibody**Synonyms:** CD203c;NPP3;PD-Ibeta**Description:** Anti-ENPP3(Xencor) Monoclonal Antibody**Background:** Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.**Applications:** Flow Cyt**Recommended Dilutions:** Flow Cyt 1:100**Host Species:** Humanized**Isotype:** IgG1**Purification:** Purified from cell culture supernatant by affinity chromatography**Species Reactivity:** Human ENPP3**Constituents:** Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.**Storage & Shipping:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).

Figure 1. Flow cytometry analysis with 1 µg/mL Anti-ENPP3(tralokinumab biosimilar) mAb (28143) on Expi293 cells transfected with Human ENPP3 protein (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram).

