

Product Description

Pioneering GTPase and Oncogene Product Development since 2010

Arl3

Catalog Number: 10152

Synonyms: ADP-ribosylation factor-like 3, ARFL3

Background: Arl3 belongs to the member of the Arf family of regulatory GTPases, within the Ras superfamily of GTPases. Arl3 is also a member of the ADP-ribosylation factor family of GTP-binding proteins. Arl3 binds guanine nucleotides but lacks ADP-ribosylation factor activity.

Amino Acid Sequence: (1-182)

MGLLSILRKLKSAPDQEVRILLLGLDNAGKTTLLK QLASEDISHITPTQGFNIKSVQSQGFKLNVWDIGG QRKIRPYWKNYFENTDILIYVIDSADRKRFEETGQ ELAELLEEEKLSCVPVLIFANKQDLLTAAPASEIAE GLNLHTIRDRVWQIQSCSALTGEGVQDGMNWV CKNVNAKKK

Source: Human, recombinant full length, His6-tag

Expression Host: E. coli Molecular Weight: 20 kDa Purity: >95 % by SDS-PAGE

Constituents: 20 mM Tris-HCl, pH 8.0, 150 mM NaCl.


Physical Appearance: White or clear

Concentration: 1 mg/mL

Storage: -80°C

Preparation Instructions

Adding of 10 mM β -mercaptoethanol or 1 mM DTT into the solution to protect the protein is recommended and using of non-ionic detergents such as n-Dodecyl β -D-maltoside (DoDM) or polyethylene detergents (e.g. C12E10) also help to stabilize the protein. Avoid repeated freezing and thawing.



The purity of His-tagged Arl3 was determined by SDS- PAGE and Coomassie Brilliant Blue Staining.

References

- Cavenagh, M. M. et al., J. Biol. Chem. 269: 18937-18942, 1994.
- 2. Grayson, C. et al., Hum. Molec. Genet. 11: 3065-3074, 2002.
- 3. Kim, H.-S. Cytogenet. Cell Genet. 83: 246-only, 1998.
- 4. Schrick, J. JAm. J. Path. 168: 1288-1298, 2006.
- 5. Veltel, S. et al., Nature Struct. Molec. Biol. 15: 373-380, 2008.