

Rab18

Catalog Number: 21050

Gene Symbol: RAB18

Description: Anti-Rab18 Rabbit Polyclonal Antibody

Background: Rab18 is a member of the Ras-related small GTPases. Rab18 plays a role in apical endocytosis/recycling.

Immunogen: A synthetic peptide from the C-terminal end of human Rab18 protein

Tested Applications: ELISA, WB, IHC

Recommended Dilutions:

ELISA 1:100-1:5000

WB 1:100-1:1000

IHC 1:100-1:1000

Concentration: 0.2 mg/ml

Host: Rabbit

Clonality: Polyclonal

Purity: Purified from serum

Format: Liquid


Preservative: No

Constituents: PBS (without Mg^{2+} and Ca^{2+}), pH7.4, 150 mM NaCl, 50% glycerol

Species Reactivity: Rab18 antibody recognizes Rab18 proteins from vertebrates.

Storage Conditions: Store at $-20^{\circ}C$. Avoid freeze / thaw cycles

Western blot:



WB: anti-Rab18 pAb

EC-1 (lane 1), HepG (lane 2) and HEK-293T (lane 3) whole-cell lysates were used for Western blots with anti-Rab18 polyclonal antibody (Cat. # 21050).

FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC APPLICATIONS