

## Rab8A

**Catalog Number:** 21071

**Gene Symbol:** RAB8A

**Description:** Anti-Rab8A Rabbit Polyclonal Antibody

**Background:** Rab8A is involved in polarized trafficking of vesicles from the *trans*-Golgi network and is the only RAB-family protein identified in mammalian cilia.

**Immunogen:** A synthetic peptide from the N-terminal region of human Rab8A

**Tested Applications:** ELISA, WB, IHC

**Recommended Dilutions:**

ELISA	1:100-1:5000
WB	1:100-1:1000
IHC	1:100-1:1000

**Concentration:** 0.2 mg/ml

**Host:** Rabbit

**Clonality:** Polyclonal

**Purity:** Purified from serum

**Format:** Liquid


**Preservative:** No

**Constituents:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH7.4, 150 mM NaCl, 50% glycerol

**Species Reactivity:** Rab8A antibody recognizes Rab8A from vertebrates.

**Storage Conditions:** Store at -20°C. Avoid freeze / thaw cycles

### Western blot:



HepG cell lysate was used for Western blot with anti-Rab8A polyclonal antibody (Cat. # 21071).

FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC APPLICATIONS