

Product Description

Pioneering GTPase and Oncogene Product Development since 2010

ERBB2

Catalog Number: 26080

Gene Symbol: ERBB2/HER2/Neu

Description: Anti-ERBB2 Rabbit Polyclonal Antibody

Background: The ErbB family is composed of four plasma membrane-bound receptor tyrosine kinases. Amplification or over-expression of the ERBB2 gene occurs in approximately 30% of breast cancers. HER2 is encoded by ERBB2, a known proto-oncogene located at the long arm of human chromosome 17(17q21-q22). HER2 proteins have been shown to form clusters in cell membranes that may play a role in tumorigenesis.

Immunogen: A synthetic peptide from the central of

ErbB2 protein.

Tested Applications: ELISA, WB, IF, IHC

Recommended Dilutions:

ELISA: 1:500-1:5000 WB: 1:200-1:2000 IHC: 1:20-1:200

Concentration: 0.2 mg/ml

Host: Rabbit

Clonality: Polyclonal

Purity: Purified from ascites

Format: Liquid
Preservative: No

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH7.4, 150

mM NaCl, 50% glycerol


Species Reactivity: Recognizes ERBB2 from

vertebrates.

Storage Conditions: Store at -20°C. Avoid freeze / thaw

cycles.

Western blot:



WB: Anti-ERBB2 pAb

Western blot analysis of ErbB2 protein.


Purified full length His-tagged ErbB2 protein (Cat. # 10208) was blotted with anti-ErbB2 rabbit polyclonal antibody (Cat. # 21202).

Web:



Pioneering GTPase and Oncogene Product Development since 2010

Immunohistochemistry



Immunohistochemical analysis of paraffin-embedded human breast cancer tissue with anti-ERBB2 polyclonal antibody (Cat. # 21202).

Tissue samples were fixed with formaldehyde and blocked with 1% serum for 15 min at 37 °C. Antigen retrieval was by heat mediation in citrate buffer (pH6). Samples were then incubated with primary antibody (1:50 dilution) overnight at 4°C. A HRP-conjugated Goat anti-Rabbit IgG (1:50) was used as secondary antibody.

Web: