

RAB3A

Catalog Number: 26031

Gene Symbol: RAB3A

Description: Anti-Rab3A Mouse Monoclonal Antibody

Background: Rab3A is a member of the Rab-family GTPases and is specifically targeted to synaptic vesicles.

Immunogen: Recombinant full length protein of Rab3A.

Tested Applications: ELISA, WB, IHC, IP

Recommended Dilutions:

ELISA	1:100-1:5000
WB	1:100-1:2000
IHC	1:100-1:1000

Concentration: 1 mg/ml

Host: Mouse

Clonality: Monoclonal

Isotype: IgG

Purity: Purified from ascites

Format: Liquid


Preservative: No

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH7.4, 150 mM NaCl, 50% glycerol

Species Reactivity: Rab3A antibody recognizes Rab3A of vertebrates.

Storage Conditions: Store at -20°C. Avoid freeze / thaw cycles

Western blot:



Wb: anti-Rab3A mAb

Western blot analysis of recombinant GST-Rab3A proteins.

FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC APPLICATIONS