

Product Description

Pioneering GTPase and Oncogene Product Development since 2010

EYA2

Catalog Number: 26096 Gene Symbol: EYA2

Description: Anti-EYA2 Mouse Monoclonal Antibody

Background: EYA2 may be involved in development of the eye. It promotes efficient DNA repair by dephosphorylating H2AX, promoting the recruitment of DNA repair complexes containing MDC1. Its function as histone phosphatase probably explains its role intranscription regulation during organogenesis.

EYA2 coactivates SIX1 and seems to coactivate SIX2, SIX4 and SIX5. It seems to be involved in myogenesis together with SIX1 and DACH2.

Immunogen: Human EYA2 protein **Tested Applications:** ELISA, WB, IHC

Recommended Dilutions:

ELISA: 1:1000-1:5000 WB: 1:100-1:1000 IHC: 1:50-1:100

Concentration: 0.8 mg/ml **Host:** Mouse

Clonality: Monoclonal

Purity: Purified from ascites

Format: Liquid Preservative: No

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH

7.4, 150 mM NaCl, 50% glycerol


Species Reactivity: Anti-EYA2 antibody

recognizes EYA2 of vertebrates.

Storage Conditions: Store at -20°C. Avoid

freeze / thaw cycles

Western blot:



WB: Anti-EYA2 mAb

Western blot analysis of recombinant full length EYA2 protein.

Purified His-tagged EYA2 protein (Cat. # 10195) was blotted with anti-EYA2 monoclonal antibody (Cat. # 26096).


Web:



Product Description

Pioneering GTPase and Oncogene Product Development since 2010

Immunohistochemistry:



Immunohistochemical analysis of paraffin-embedded thyroid cancer tissue with anti-EYA2 monoclonal antibody (Cat. #26096).

Tissue samples were fixed with formaldehyde and blocked with 1% serum for 15 min at 37 °C. Antigen retrieval was by heat mediation in citrate buffer (pH6). Samples were then incubated with primary antibody (1:50) overnight at 4°C. A HRP-conjugated Goat anti-mouse IgG (dilution 1:50) was used as secondary antibody.

Web: