

Product Description

Pioneering GTPase and Oncogene Product Development since 2010

EGFR(G719A)

Catalog Number: 26115

Gene Symbol: EGFR/HER1/ERBB1

Description: Anti-EGFR (G719A) Mouse Monoclonal Antibody

Background: EGFR is a receptor tyrosine kinase binding

ligands of the EGF family and activating several signaling cascades to convertextracellular cues into appropriate cellular responses. It activates at least 4 major downstream signaling cascades including the RAS-RAF-MEK-ERK, PI3 kinase-AKT, PLCgamma-PKC and

STATs modules. Defects in EGFR are associated with lung cancer (LNCR).

Immunogen: A synthetic peptide from the internal region of EGFR which includes themutation of G719A, human

origin.

Tested Applications: ELISA, WB, IF, IHC

Recommended Dilutions:

ELISA: 1:1000-1:5000 WB: 1:500-1:1000 IHC: 1:50-1:100

Concentration: 1 mg/ml

Host: Mouse

Clonality: Monoclonal

Isotype: IgG

Purity: Purified from ascites

Format: Liquid Preservative: No

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH7.4, 150

mM NaCl, 50% glycerol


Species Reactivity: Recognizes EGFR (G719A), but not

wild type EGFR from vertebrates.

Storage Conditions: Store at -20°C. Avoid repeated

freezing and thawing

Western blot:



WB: Anti-EGFR (G719A) mAb

Western blot analysis of recombinant EGFR(G719A) and wild type proteins.

Purified His-tagged EGFR (G719A) protein (aa 702-987) and wild type protein (aa 702-987) were blotted with anti-EGFR(G719A) mouse antibody (Cat. #


26115).

Web:



Pioneering GTPase and Oncogene Product Development since 2010

Immunofluorescence:



Immunofluorescence of cells expressing EGFR proteins with anti-EGFR (G719A) antibody.

HEK293T cells were transfected with pCDNA3-GFP-EGFR (WT) plasmid (left column) or pCDNA3-GFP-EGFR (G719A) plasmid (right column), then fixed and stained with anti- EGFR (G719A) monoclonal antibody (Cat. # 26115).